


Spare Parts List Babytherm 8000 WB/OC, US version




Summary of pictures:

Picture	Description
1	BABYTHERM 8000 OC
2	BABYTHERM 8000 WB
3	TROLLEY, FOR HEIGHT ADJUSTMENT
4	PEDAL FOR HEIGHT ADJUSTMENT
5	PLATE FOR HEIGHT ADJUSTMENT (CABLE HARNESS)
6	TROLLEY WITH FIXED HEIGHT BABYTHERM 8000 OC/WB
7	BASE PLATE, LARGE, WITH HANDLE
8	HANDLE
9	BASE PLATE, SMALL, WITH HANDLE
10	HANDLE
11	ELECTRIC OUTLET PANEL WITH POWER CORD
12	INSTALLATION SET CUPBOARD BABYTHERM
13	SWIVEL CUPBOARD, COMPLETE
14	TRIPOD WITH TRAY
15	DEPOSIT PLATE
16	HANDLE, COMPLETE / FOR DEPOSIT PLATE
17	COT BABYTHERM 8000
18	INCLINATION, COMPLETE
19	LOCK, COMPLETE.
20	SIDE WALLS, HEIGHT 15CM
21	SIDE WALLS, HEIGHT 23CM
22	HOOD BABYTHERM 8000
24	TABLE INCLINABLE, COMPLETE
25	MATTRESS HEATING FOR GEL MATTRESS
26	CONTROL UNIT, COMPLETE, PART 1
27	CONTROL UNIT, COMPLETE, PART 2
28	MODIFICATION KIT 100/127V
29	HEATING PLATE, COMPLETE
30	TEMPERATURE DETECTOR, COMPLETE


BABYTHERM 8000 OC

Picture 1



Item No.	Description	Order-Code
2	MATTRESS (FOAM, FOR UNITS W/O MATTRESS HEATING ONLY)	2M21012
2	GEL-MATTRESS (FOR UNITS WITH MATTRESS HEATING ONLY)	2M20827
	REPAIR-KIT FOR GEL-MATTRESS	2M21324

Item No.	Description	Order-Code
3	INSTALLATION SET CUPBOARD, SEE PICTURE 12	
4	SWIVEL CUPBOARD, SEE PICTURE 13	
5	PEDAL, SEE PICTURE 4	
6	TROLLEY, LARGE FOR HEIGHT ADJUSTMENT, SEE PICTURE 3	
6	TROLLEY, SMALL FOR HEIGHT ADJUSTMENT, SEE PICTURE 3	
7	ELECTRIC OUTLET PANEL 110/127V, WITH POWER CORD, SEE PICTURE 11	
8	ALTITUDE ADJUSTMENT, 110V / 86-116CM, WHEN REPLACING HEIGHT ADJUSTMENT ORDER IN ADDITION LABEL SET USA P/N 2M21896	2M20946
8	ALTITUDE ADJUSTMENT, 110V / 92-122CM, WHEN REPLACING HEIGHT ADJUSTMENT ORDER IN ADDITION LABEL SET USA P/N 2M21896	2M21175
9	BASE PLATE, LARGE WITH HANDLE, SEE PICTURE 7	
10	MATTRESS HEATING, GEL-MATTRESS, SEE PICTURE 25	
12	SIDE WALLS, HEIGHT 15CM, SEE PICTURE 20	
12	SIDE WALLS, HEIGHT 23CM, SEE PICTURE 21	
13	TRIPOD WITH TRAY	2M21066
17	INTERIOR DISC	2M20936
18	FIXED HEIGHT TROLLEY, BT 8000 OC, SEE PICTURE 6	
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896


BABYTHERM 8000 WB**Picture 2**

Item No.	Description	Order-Code
1	SIDE WALLS, HEIGHT 15CM, SEE PICTURE 20	
1	SIDE WALLS, HEIGHT 23CM, SEE PICTURE 21	
2	INTERIOR DISC	2M20936

Item No.	Description	Order-Code
4	MATTRESS (FOAM, FOR UNITS W/O MATTRESS HEATING ONLY)	2M21012
4	GEL-MATTRESS (FOR UNITS WITH MATTRESS HEATING)	2M20827
	RAPAIR-KIT FOR GEL-MATTRESS	2M21324
5	INSTALLATION SET CUPBOARD, SEE PICTURE 12	
6	ELECTRIC OUTLET PANEL, WITH POWER CORD, 110/127V, SEE PICTURE 11	
7	SWIVEL CUPBOARD, SEE PICTURE 13	
8	PEDAL, SEE PICTURE 4	2M20290
9	TROLLEY, SMALL FOR HEIGHT ADJUSTMENT, SEE PICTURE 3	
9	TROLLEY, LARGE FOR HEIGHT ADJUSTMENT, SEE PICTURE 3	
10	ALTITUDE ADJUSTMENT, 110V / 86-116CM, WHEN REPLACING HEIGHT ADJUSTMENT ORDER IN ADDITION LABEL SET USA P/N 2M21896	2M20946
10	ALTITUDE ADJUSTMENT, 110V / 92-122CM, WHEN REPLACING HEIGHT ADJUSTMENT ORDER IN ADDITION LABEL SET USA P/N 2M21896	2M21175
11	BASE PLATE, SMALL WITH HANDLE, SEE PICTURE 9	
12-13	CONTROL UNIT, SEE PICTURE 26	
14	SUPPORT	2M21342
15	HOOD BABYTHERM 8000, SEE PICTURE 22	2M21030
16	TROLLEY, FIXED HEIGHT, BT 8000 WB, SEE PICTURE 6	
16	TROLLEY, FIXED HEIGHT PEDESTAL, BT 8000 WT, (without illustration), SEE PICTURE 6	
17	VENTILATION HOSE SUPPORT, (without illustration)	2M21191
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896


TROLLEY, FOR HEIGHT ADJUSTMENT

Picture 3




Item No.	Description	Order-Code
	TROLLEY, LARGE, FOR HEIGHT ADJUSTMENT, WHEN REPLACING HEIGHT ADJUSTEMENT ORDER IN ADDITION LABEL SET USA P/N 2M21896	2M21036
1	TUBE, RIGHT	2M20951
2	TUBE, LEFT	2M20950
	TROLLEY, SMALL, FOR HEIGHT ADJUSTMENT, WHEN REPLACING HEIGHT ADJUSTEMENT ORDER IN ADDITION LABEL SET USA P/N 2M21896	
5	CAP	2M20348
6	FIXING GUIDE ROLL 125 MM	2M20794
7	GUIDE ROLL 125 MM	2M20792
8	CH. HEAD SCREW M10X40 DIN912-8.8	1308564
12	CYL. HEAD SCREW M 6X25 DIN 912	1325779
13	LOCK WASHER	G12768
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896


Picture 4



Item No.	Description	Order-Code
1-8	PEDAL, COMPLETE.	2M20290
1	PEDAL	2M20218
4	PLATE (CABLE HARNESS), SEE PICTURE 5	

PLATE FOR HEIGHT ADJUSTMENT (CABLE HARNESS)
Picture 5



Item No.	Description	Order-Code
1-7	PLATE, CPL.	2M20370
2	MICROSWITCH SDRU 1POL TAWE	1831666
5	FLAT PLUG DIN 46244-A6,3	1814079
6	SPLIT WASHER 4 DIN7980-X12CRNI	1293680
7	SCREW AM 4X10 DIN84-A4/051	1329049

TROLLEY WITH FIXED HEIGHT BABYTHERM 8000 OC/WB**Picture 6**


Item No.	Description	Order-Code
1 - 5	TROLLEY FIXED HEIGHT BT 8000 OC, LOW, WHEN REPLACING COLUMN ORDER IN ADDITION LABEL SET USA P/N 2M21896	2M22018
3	CAP	2M20989
4	FIXING CASTER	2M21050
5	CASTER	2M21048
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896

BASE PLATE, LARGE, WITH HANDLE


Picture 7




Item No.	Description	Order-Code
1-7	BASE PLATE, LARGE, WITH HANDLE, WHEN REPLACING BASE PLATE ORDER IN ADDITION LABEL SET USA BABYTHERM 8000/8004/8010 P/N 2M21896	2M21146
2	HANDLE, SEE PICTURE 8	
4	SCREW M6X6 DIN 916-A4	1320831
5	FLAT PLUG DIN 46244-A6,3	1814079
6	CH. HEAD SCREW M4X8 DIN912	1263056
7	SPLIT WASHER B4 DIN127-X12CRNI	1331345
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896

HANDLE**Picture 8**


Item No.	Description	Order-Code
2	END PIECE	2M20696
3	TUBE	2M13388
4	CAP SCREW M 8X75 DIN 912	1326473

BASE PLATE, SMALL, WITH HANDLE
Picture 9



Item No.	Description	Order-Code
1-5	BASE PLATE, SMALL, WITH HANDLE, WHEN REPLACING BASE PLATE ORDER IN ADDITION LABEL SET USA BABYTHERM 8000/8004/8010 P/N 2M21896	2M21148
2	HANDLE, SEE PICTURE 10	
3	FLAT PLUG DIN 46244-A6,3	1814079
4	CH.HEAD SCREW M4X8 DIN912	1263056
5	SPLIT WASHER B4 DIN127-X12CRNI	1331345
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896

HANDLE**Picture 10**


Item No.	Description	Order-Code
2	CLOSING PLUG	G60036
3	TUBE	2M13388
4	CAP SCREW M 8X75 DIN 912	1326473

INSTALLATION SET CUPBOARD BABYTHERM**Picture 12**


Item No.	Description	Order-Code
1-4	INSTALLATION SET CUPBOARD BABYTHERM	2M20868
3	SCREW DIN 912-M8X390+-2 8.8/155	1342096

SWIVEL CUPBOARD, COMPLETE
Picture 13



Item No.	Description	Order-Code
1-6	SWIVEL CUPBOARD, CPL.	2M20638
1	SWIVELLING DRAWER	2M20627
2	INSERT	2M20642
3	SWIVELLING PLATE	2M20628
4	TUBE	2M20639
5	CATCH	2M20632
6	SPRING	2M20645

TRIPOD WITH TRAY**Picture 14**

Item No.	Description	Order-Code
2	CABLE CLIP (3,6X140)	8300358
4	SCREW M6X6 DIN 916-A4	1320831
5	HEXAGON NUT M4 DIN 934-A4/051	1328956
6	SPLIT WASHER B4 DIN127-X12CRNI	1331345
7	SCREW AM 4X10 DIN84-A4/051	1329049
9	GASKET	2M19591
11	DEPOSIT PLATE, SEE PICTURE 15	

DEPOSIT PLATE**Picture 15**


Item No.	Description	Order-Code
1-6	DEPOSIT PLATE	2M21158
3	SUPPORT II	2M21214
4	SUPPORT I	2M21160
5	BUFFER	6800478
	HANDLE, (without illustration, not included in P/N 2M21158), SEE PICTURE 16	

HANDLE, COMPLETE / FOR DEPOSIT PLATE**Picture 16**


Item No.	Description	Order-Code
1-3	HANDLE, CPL.	2M21161
2	END PIECE	2M20696
3	CAP	G13569

COT BABYTHERM 8000

Picture 17




Item No.	Description	Order-Code
1	COT (WITH CORNERS), WHEN REPLACING COT ORDER IN ADDITION LABEL SET USA P/N 2M21896	2M20832
2	LOCK, COMPLETE, SEE PICTURE 19	
3	INCLINATION, SEE PICTURE 18	
6	CORNER WITH SURFACE FINISH	2M21390
	SET HOOKED RAIL (without illustration)	2M21150
	HOOK (without illustration)	2M21293
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896

INCLINATION, COMPLETE**Picture 18**


Item No.	Description	Order-Code
4	TOROIDAL SEALING RING	M12747

LOCK, COMPLETE

Picture 19




Ansicht Z



Item No.	Description	Order-Code
1 - 16	LOCK, CPL.	2M20838
3	STOP BUFFER	2M13141
4	LOCK WASHER	G12768
8	CABLE BUSHING	CH19989
16	CATCH	2M20855

SIDE WALLS, HEIGHT 15CM


Picture 20




Item No.	Description	Order-Code
1 - 6	SIDE WALLS, HEIGHT 15CM, ORANGE, ORDER IN ADDITION LABEL SET USA (BABYTHERM 8000/8004/8010) P/N 2M21896	2M21032
1,2,6-8	SIDE WALLS, HEIGHT 15CM, WHITE ORDER IN ADDITION LABEL SET USA (BABYTHERM 8000/8004/8010) P/N 2M21896	2M30242
1	SIDE PANE (15 CM), ORANGE, ORDER IN ADDITION LABEL SET USA (BABYTHERM 8000/8004/8010) P/N 2M21896	2M20886
2	FRONT PANE (15CM)	2M20891
3	REAR PANE (15 CM)	2M20892
4	REPAIR SET HINGE COMPLETE, RIGHT, ORANGE	2M21082
5	REPAIR SET HINGE COMPLETE, LEFT, ORANGE	2M21084
6	TUBE BUSHING	2M20434
7	SPRING-LOADED CATCH, white	2M30210
7	SPRING LOADED CATCH, orange	2M20898
8	BUTTON, white	2M30211
8	BUTTON, orange	2M20899
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896

SIDE WALLS, HEIGHT 23CM


Picture 21




Item No.	Description	Order-Code
1 - 8	SIDE WALLS, HEIGHT 23CM, ORDER IN ADDITION LABEL SET USA (BABYTHERM 8000/8004/8010) P/N 2M21896	2M21034
1-3, 6-10	SIDE WALLS, HEIGHT 23CM, WHITE ORDER IN ADDITION LABEL SET USA (BABYTHERM 8000/8004/8010) P/N 2M21896	2M30301
1	SIDE PANE (23 CM), ORANGE, ORDER IN ADDITION LABEL SET USA (BABYTHERM 8000/8004/8010) P/N 2M21896	2M20931
2	FRONT PANE (23 CM), ORANGE	2M20932
3	REAR PANE (23 CM), WITH 3 OPENINGS, ORANGE	2M20937
4	REPAIR SET HINGE, RIGHT, ORANGE	2M21082
5	REPAIR SET HINGE, LEFT, ORANGE	2M21084
6	TUBE BUSHING	2M20434
7	LOCK	2M13042
8	LABEL	2M21147
9	SPRING-LOADED CATCH, white	2M30210
9	SPRING LOADED CATCH, orange	2M20898
10	BUTTON, white	2M30211
10	BUTTON, orange	2M20899
	LABEL SET USA BABYTHERM 8000/8004/8010	2M21896

HOOD BABYTHERM 8000**Picture 22**

Item No.	Description	Order-Code
1 - 8	HOOD BABYTHERM 8000	2M21030
6	THERMOMETER	2M11111
7+8	SET OF SUPPORTS FOR THERMOMETER	2M16752

TABLE INCLINABLE, COMPLETE
Picture 24



Item No.	Description	Order-Code
1 - 8	TABLE INCLINABLE, COMPLETE	2M21186
4	SCREW DWN 561-4X20-C-A2-H	1338072
7	WASHER A 8,4 DIN 125-A4/051	1324780
9 - 11	CLAW, COMPLETE	2M21378
10	THRUST WASHER 5 DIN 6799	1326406
11	ADJUSTING SCREW	M20647

MATTRESS HEATING FOR GEL MATTRESS**Picture 25**

Item No.	Description	Order-Code
1	CONTROL UNIT, SEE PICTURE 26	
3	GEL-MATTRESS / BABYTHERM 8000	2M20827

CONTROL UNIT, COMPLETE, PART 1


Picture 26




Item No.	Description	Order-Code
	SERVICE EXCHANGE CONTROL UNIT (2M21500)	SE2M21506
	MODIFICATION KIT 100/127V (FOR 2M21500), SEE PICTURE 28	
1	FUSIBLE PLUG T1,6A DIN41662 (100/127V UNITS)	1815105
2	CABLE HARNESS MAINS	8201795
4	BUTTON, BLACK	8200619
5	KEY PAD, CPL.	2M20885
6	PCB WT-CONTROLLER-MATTRESS HEATING.	8290668
7A	PCB POWER (OLDER VERSION, REFER TO TSB BABYTHERM 8000 # 4)	8290641
7B	PCB POWER (NEWER VERSION, REFER TO TSB BABYTHERM 8000 # 4)	8290771
10	KREL 5V/8A	1837001
11	FUSE LINK DIN 41662 T6,3A	1815172
14	FLAT PLUG DIN 46244-A6,3	1814079
15	XDSUB LOCK BOLT	1822357
16	PCB DISPLAY	8290611
	D PROG. KEY	1820877
17	SOFTWARE BABYTHERM 8000, EPROM 1.04	2M22467
18	FERRITE BUSHING DM6-7	1834967

CONTROL UNIT, COMPLETE, PART 2

Picture 27




Item No.	Description	Order-Code
2	BUTTON, BLACK	8200619
4	HOSE 7X2,5 SI NF M29908	1198343
6	ISOLATION TRANSFORMER	8201796

MODIFICATION KIT 100/127V
Picture 28



Item No.	Description	Order-Code
1 - 4	MODIFICATION KIT 100/127V	2M21176
2	FUSIBLE PLUG T1,6A DIN41662	1815105

HEATING PLATE, COMPLETE

Picture 29



Item No.	Description	Order-Code
1-12	HEATING PLATE, COMPLETE	2M21510
3	PROFILE PACKING	2M21423
5	CABLE CLIP (2,4X92)	8712007
6	SCREW AM 4X5 DIN 85-A4	1330063
7	SPLIT WASHER B4 DIN127-X12CRNI	1331345
8	TEMPERATURE DETECTOR, SEE PICTURE 30	
12	RECESSED HEAD SCREW M4X20 DIN7985	1341308

TEMPERATURE DETECTOR, 2-FOLD**Picture 30**

Item No.	Description	Order-Code
1-11	TEMPERATURE DETECTOR, COMPLETE	2M20859